Medical Gas Flowmeters

Instructions for Use

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Rx Only

Safety Instructions

This manual provides you with important information about Medical Gas Flowmeters and should be read carefully to ensure the safe and proper use of this product. READ AND UNDERSTAND ALL THE SAFETY AND OPERATING INSTRUCTIONS CONTAINED IN THIS BOOKLET BEFORE USING THIS PRODUCT. IF YOU DO NOT UNDERSTAND THESE INSTRUCTIONS, OR HAVE ANY QUESTIONS, CONTACT YOUR SUPERVISOR, DEALER OR THE MANUFACTURER BEFORE ATTEMPTING TO USE THE DEVICE.

Intended Use

A compensated thorpe tube flowmeter is a device intended for medical purposes that is used to control and measure gas flow rate accurately. The device includes a vertically mounted tube with the outlet of the flowmeter calibrated to a reference pressure.

WARNING	Indicates a potentially hazardous situation, if not avoided, could result in death or serious injury.				
ATTENTION	Indicates a potentially hazardous situation, if not avoided, could result in minor or moderate injury.				
CAUTION	Indicates a potentially hazardous situation, which, if not avoided could result in damage to the device or other property.				
	Caution, consult accompanying documents			Read L/min at center of ball	
	Use no oil		Manufacturer		Serial Number
	Manufacturer's authorized representative in the European Union				
	Catalog Number (Device Identifier). This identifier includes alpha-numeric characters that correspond to the vacuum regulator model, gauge type, color and any fittings ordered. It is located on the outer package label of your unit.				
	Symbol indicates the device complies with the requirements of Directive 93/42/EEC concerning medical devices (on CE marked devices only)				
Rx Only	Caution: Federal (USA) Law restricts this device to sale by or on the order of a licensed healthcare provider.				

Receiving/Inspection

Remove product from package and inspect for damage. If product is damaged, DO NOT USE and contact your dealer or equipment provider.

CAUTION The tube inside the 0-200 cc/min and 0-1 L/min flowmeter is made of glass and is fragile. Special care should be taken to avoid breaking the glass.

ATTENTION It is very important to allow product to remain in original packaging for 24 hours to acclimate to room temperature before use.

ATTENTION Store the product in a sealed package to avoid environmental damage. The operating and storage temperature for the Flowmeter should reflect typical environmental conditions of a medical facility environment.

User Responsibility

WARNING Service of this device should ONLY be performed by properly trained individuals. The Flowmeter must ONLY be used under the direction of a healthcare professional who is trained in its use.

This product performs as explained in this manual, as long as the assembly, use, repair and maintenance are properly followed according to our instructions. Periodic review of this device is recommended. If any damage or defects are present, the product should not be used. This includes parts that may have been altered, contaminated, worn or missing. If any of the above are noted, immediate repair / replacement is required. In compliance with the manufacturers warranty, repair of this device is not to be performed by anyone other than a qualified professional. If this device is subject to improper maintenance, repair, use and/or abuse leading to malfunction of the device, replacement is the sole responsibility of the user.

WARNING Due to the possibility of explosion caused by static charge, operation of this device is not to be done if flammable anesthetics are present.

TO MINIMIZE THE RISK OF EXPLOSION OR FIRE:

- NEVER attempt to attach a Flowmeter directly to a cylinder
- NEVER use grease, oil, organic lubricants or flammable materials on or near the Flowmeter
- NEVER smoke in an area where oxygen is being used
- NEVER use any type of flame or flammable or explosive material near the Flowmeter
- ALWAYS follow CGA and ANSI standards for Flowmeters and Medical Gas Products (E-7) and Oxygen Handling (G-4)

ATTENTION Keep the Flowmeter in a clean area when not being used.

ATTENTION Ensure that all the connections are tightened and free of leaks prior to use. Use only an oxygen-safe leak detector when testing for leaks.

WARNING Each Flowmeter is for use with only one type of gas.

Specifications*

Scale	Increments	Accuracy Averages	Max. Flood/Flush	Transport/Storage Requirements
0-200 cc/min**	25 cc/min (starts at 25 cc/min)	+/- 20 cc/min	500 cc/min	-40°F to 140°F (-40°C to 60°C)
0-1 L/min**	0.1 L/min (starts at 0.1 L/min)	+/- 0.1 L/min	2.5 L/min	-40°F to 140°F (-40°C to 60°C)
0-3.5 L/min**	0.125 L/min (from 0.125-1 L/min) 0.25 L/min (from 1-3.5 L/min)	+/- 0.5 L/min above 0.5 L/min	45 L/min	-40°F to 140°F (-40°C to 60°C)
Scale	Increments	Accuracy Averages	Min. Flood/Flush	Transport/Storage Requirements
0-8 L/min	0.5 L/min (starts at 0.5 L/min)	+/- 0.5 L/min or +/- 10% of reading (whichever is greater)	50 L/min	-40°F to 140°F (-40°C to 60°C)
0-12 L/min	0.5 L/min (from 1-6 L/min) 1 L/min (from 6-12 L/min)	+/- 0.5 L/min or +/- 10% of reading (whichever is greater)	50 L/min	-40°F to 140°F (-40°C to 60°C)
0-15 L/min	0.5 L/min (from 1-5 L/min) 1 L/min (from 5-15 L/min)	+/- 0.5 L/min or +/- 10% of reading (whichever is greater)	50 L/min	-40°F to 140°F (-40°C to 60°C)
0-16 L/min	0.5 L/min (from 1-6 L/min) 1 L/min (from 6-16 L/min)	+/- 0.5 L/min or +/- 10% of reading (whichever is greater)	50 L/min	-40°F to 140°F (-40°C to 60°C)
0-30 L/min	2 L/min (from 4-30 L/min) (starts at 3 L/min)	+/- 0.5 L/min or +/- 10% of reading (whichever is greater)	70 L/min	-40°F to 140°F (-40°C to 60°C)
0-70 L/min	5 L/min (starts at 10 L/min)	+/- 10% of reading	75 L/min	-40°F to 140°F (-40°C to 60°C)

*Specifications are nominal, subject to change without notice. Flowmeters are calibrated at the pressure indicated on the flow tube, 70°F (21°C), at standard atmospheric pressure. **For use in Neonatal, Pediatric and other low flow applications.

WARNING Prior to administering medical air or oxygen, confirm required flow and frequently monitor the flow.

WARNING 0-200 cc/min, 0-1 L/min, and 0-3.5 L/min Flowmeters are NOT to be used for resuscitation.

MRI WARNING This product may contain magnetic, ferrous material that may affect the result of an MRI.
MR conditional options may be available. Contact your local Ohio Medical Sales Representative.

Operating Instructions

1. Turn Flowmeter off by turning knob fully clockwise.
2. Inspect the Flowmeter for damage. If any damage is found, DO NOT use the Flowmeter.

CAUTION Over tightening the knob when turning the Flowmeter off will cause damage. The Flowmeter **MUST** be positioned vertically to ensure maximum accuracy.

3. Connect the Flowmeter to the supply pressure specified on the Flow Tube.

NOTE: Supply pressure of 50 psi or 60 psi.

WARNING The accuracy may be affected if the temperature of the gas is different than 70°F (21°C) and the supply pressure is different than indicated on the Flow Tube.

NOTE: The accuracy of the flow will not be affected by the attachment of accessories, however, the indicated flow may change.

WARNING Connection to the gas source must be done by using only the appropriately indexed fitting.

4. Ensure that the Float Ball is at the very bottom of the Flow Tube when turned off.

NOTE: If the Float Ball is not at the bottom of the Flow Tube, the Flowmeter could be leaking. Please contact your dealer or Ohio Medical, LLC.

5. Adjusting the Flow: To DECREASE Flow: Turn the knob clockwise.
To INCREASE Flow: Turn the knob counter clockwise.

6. To set the flow, align the center of the Float Ball to the indicator line on the Flow Tube.

WARNING To avoid injury **ALWAYS** confirm flow requirement for patient prior to dispensing. Check flow frequently while being administered to patient.

7. If the flow is adjusted beyond the highest calibrated indicator, an undetermined flow will arise.

8. Turn knob completely counter clockwise to achieve maximum Flow/Flush Flow.

NOTE: Any flow beyond the highest calibrated line of the Flow Tube with unrestricted flow is Flood/Flush Flow.

Cleaning Instructions

Use a clean damp cloth with a mild cleaning solution to wipe outside of product. **DO NOT** gas sterilize with Ethylene Oxide (ETO). **DO NOT** clean with pungent hydrocarbons.

CAUTION **DO NOT** submerge product in any form of liquid. This will cause damage and void the warranty.

Troubleshooting

If the Flowmeter does not function properly, contact your dealer or equipment provider.

Maintenance Prevention

Before and after each use, **INSPECT THE FLOWMETER FOR DAMAGE.**

WARNING When changing connectors on the Flowmeter for service or replacement, **NEVER** re-attach connectors of a different gas. Doing so may result in patient injury or damage to the equipment

CAUTION Disconnect Flowmeter from gas supply **BEFORE SERVICING.**

WARNING Ensure shroud is securely tightened to prevent disconnection under pressure. Turn Clockwise to tighten.

Warranty

The standard warranty period is 60 months from date of receipt by the customer (with the exception of the 0-200 cc/min, 0-1 L/min, and 0-70 L/min flowmeter for 12 months). The standard warranty is only valid for medical devices handled according to the instructions for use (IFU) and general industry good practice and standards.

EC REP

EMERGO EUROPE

Westervoortsedijk 60 6827 AT Arnhem The Netherlands

AUSTRALIAN SPONSOR: EMERGO AUSTRALIA

Level 20, Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia

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