

## **Product Data Sheet**

20 Cu. Ft Standard Series Laboratory Manual Defrost Freezer

# CMEB-FRZ-20-S-HCF





### **Product Description**

Laboratory freezers offer professional cold storage for scientific and industrial settings. Units come equipped with temperature alarms, keyed door locks, forced draft circulation, probe access, and precise temperature monitoring through microprocessor control.

Our Standard Series freezers provide multi-functional features in a costeffective design. All units run on natural, hydrocarbon refrigerant for environmental health and energy efficiency.

## Certifications





**General Description and Application** Storage capacity (cu. ft) 20 Right Hinged Shelves 6 Fixed Non-applicable Baskets Non-applicable Mounting Leveling Legs Interior lighting Non-applicable Airflow management Non-applicable External probe access Rear wall port (3/4") dia. Insulation Cabinet is foamed-in-place with EPA compliant high density urethane foam Exterior materials White powder coated steel Access Control Keyed door lock General warranty Two (2) year parts and labor warranty Products require 4" of clearance around back and sides of the unit for proper ventilation - improper installation will void the warranty and will lead to temperature maintenance issues, overheating and possible failure Warranty Disclaimer Product Weight (lbs) 225 250 Shipping Weight (lbs) Rated Amperage Power Plug/Power Cord NEMA 5-15 Plug, conforms to UL471 requirements Facility Electrical Requirement 115V, 60 Hz Agency Listing and Certification ETL, C-ETL listed and certified to UL471 standard, hydrocarbon refrigerant safety, Energy Star Certified

# Refrigeration System

Refrigerant R290 Defrost Manual

# Controller, Configuration, Alarms and Monitoring

Microprocessor temperature controller with digital temperature display, °C/°F switchable Controller technology

Display technology Standard Digital Temperature Display

Adjustable Temperature Range -15°C to -25°C

State switching remote alarm contacts External alarm connection Alarms

Audible and visual high and low temperature alarms


Controller probe 2 Probes; 1 in air and 1 in sample bottle Non-applicable

Digital Communication Data Transfer Non-applicable

Defrosting of the unit will need to be performed periodically to prevent severe ice build-up

Indoor use only, +18°C to +26°C (+65°F to +78°F), <70% RH

Specifications are subject to change



Rev\_07202023